

LOCAMED CLOSED CAMERA COVER



REF.	DESCRIPTION	QTY / BOX
LM8007	Closed camera sleeve made of tough, thin polyethylene. This disposable device is sterile and for single-use only. Dim. 152 x 2400 x 0.05 mm	50 Units/box

- **Fast, simple and safe application to create a sterile field around the camera head and lead**
- **Quick, easy and secure switching of the endoscope during laparoscopic procedures**
- **Also suitable for cystoscopy or arthroscopic procedures**
- **Enhanced sterile barrier**
- **No tape required**
- **Sleeve contained in a plastic ring for easy deployment and covering of the camera head and lead**