




# Disposable Veress Needle NEW 200mm length

## Cost effective

Fairmont Medical's Veress Needle is available in three convenient sizes. The needle is spring loaded with a inner stylet and well positioned grip at the top of the needle for ease of use.

- Shaped soft touch hub, designed for better grip and handling
- Viewing window, with red indicator
- Coloured Coded 120mm, 150mm and 200mm length
- Worldwide, Fairmont Medical are the only manufacturers who sell a 200mm length Veress Needle
- The 200mm Veress Needle is geared towards bariatric surgery
- Competitively priced



Sterile

Product Code	Colour Grip	Length	Box Quantity
DVN7012	Blue	120mm	30
DVN7015	Yellow	150mm	30
DVN7200	Orange	200mm	30



DVN7012



DVN7200



DVN7015



Balmer Médical SA  
Rte de Provence 52  
CH-1426 Concise

Tél. +41 24 436 29 07  
info@balmermedical.ch  
www.balmermedical.ch

