

Surgical
Innovations


 balmer medical

The No.1 Choice for
Minimally Invasive Surgery.

INSTRUMENTATION




RETRACTION



PORT ACCESS



A faint, light gray wireframe illustration of a human figure, showing the head, neck, shoulders, and arms. The figure is composed of a grid of lines that define the body's structure, giving it a technical or medical appearance. It is positioned on the left side of the page, with the head at the top left and the arm extending downwards and to the right.

Surgical Innovations.
We listen,
We innovate,
We improve,
We care...

Quick™Range DISPOSABLE

DISPOSABLE
PRODUCT

Quick™Range

Quick™ sterile-packed disposable laparoscopic instruments are designed to deliver ultimate performance in a simple, convenient format. Each device comes complete with a handle and insert. A range of insert configurations are available.



Key Strengths

- + Disposable – maximum convenience and ultimate performance at minimum cost
- + Electro-cautery connection is positioned for unobtrusive monopolar diathermy
- + Large tactile rotation cone ensures accurate control of the distal tip
- + Ergonomically designed handle enhances tactile feedback for high-precision procedures
- + Ratchet mechanism is conveniently located for ease of actuation
- + Quick™ jaws are profiled to ensure atraumatic manipulation of tissue
- + A range of instrument configurations available – Quick™Cut (scissors) and Quick™Range (graspers and dissectors)
- + Available with familiar ‘pistol-grip’ handle with or without ratchet
- + Quick™Cut blades are precision-ground for optimum cutting performance

Quick™Range

Quick™Cut Scissor

(16mm blade option)



Product Name
Quick™Cut Scissor
Product Code
Q21201
Product Description
Without Ratchet
Blade Length
16mm

Quick™Cut Mini Scissor

(10mm blade option)



Product Name
Quick™Cut Mini Scissor
Product Code
Q21202
Product Description
Without Ratchet
Blade Length
10mm

Quick™Dissect Curved Dissector



Product Name
Quick™Dissect Curved Dissector
Product Code
Q21205
Product Description
Without Ratchet
Blade Length
15mm

Quick™Clinch Atraumatic Grasper



Product Name
Quick™Clinch Atraumatic Grasper
Product Code
Q21210
Product Description
With Ratchet
Blade Length
20mm

Quick™Grasp Fenestrated Grasper



Product Name
Quick™Grasp Fenestrated Grasper
Product Code
Q21218
Product Description
With Ratchet
Blade Length
20mm

Delivering products designed, developed and manufactured in the UK



 balmer medical

Surgical
Innovations

Balmer Médical SA
Rte de Provence
1426 Concise
info@balmermedical.ch
Tél. +41 24 436 29 07

www.balmermedical.ch

Leaflet Part Number: ALL005 Issue1
Change Request Number: CR3968