

THD N-Ano



THD N-Ano is a disposable self-illuminating diagnostic anoscope with a powerful and fully integrated LED lighting system for an optimal visualization of the anal canal.

Thanks to its small diameter, the use of this anoscope is recommended on patients presenting a sphincter hypertone or stenosis. The comfortable handle design makes it particularly easy to use. The special design helps better direct the light beam to the area being inspected

BENEFITS

- Self-illuminating with in-built lighting system
- Ready to use
- No cross contamination risks
- Lightweight
- Small diameter
- Frosted end tip for an enhanced illumination
- Comfortable handle design
- No sterilization costs
- Individually packed

THD N-Ano

Length 65 mm
Diameter 18 mm
Code Bevelled: 800134
Truncated: 800040

THD N-Ano Maxi

Length 65 mm
Diameter 23 mm
Code Bevelled: 800166
Truncated: 800165

