


Surgical
Innovations

 balmer medical

The No.1 Choice for
Minimally Invasive Surgery.

RETRACTION




INSTRUMENTATION



PORT ACCESS



A wireframe illustration of a human figure, showing the head, neck, and upper torso. The figure is composed of a grid of thin, light-colored lines that define the contours of the body. The background is plain white.


Surgical Innovations.
We listen,
We innovate,
We improve,
We care...

PretzelFlex™ REUSABLE



PretzelFlex™ 3mm and 5mm

PretzelFlex™ is the world's first reusable, pretzel-shaped organ and tissue retractor for use in Minimally Invasive Surgery (MIS). It provides improved access and visualisation of surgical sites. The device has a slimline shaft that can be inserted through a 3mm or 5mm port. Then, once inside the abdomen, the PretzelFlex™ forms into the innovative pretzel shape.




Key Strengths

- + Fully reusable
- + Available in 3mm or 5mm variations
- + Impressive strength and stability to support large organs
- + Large surface area ensures increased contact which reduces pressure on the retracted organs
- + Unique pretzel shape provides better access and visualisation of surgical sites
- + Creates shape in smaller space envelope compared with triangular retraction
- + Easily inserted through either YelloPort+plus™ 3mm or YelloPort+plus™ 5mm Resposable® port access system

'Advance and Actuate' Technique

PretzelFlex™ uses SmartFlex™ technology, which allows the actuation of the device to begin before all the segments have entered the body. This 'easy-to-use' technique minimises the space required to form the innovative pretzel shape and means the shape can be created within a smaller space envelope than competitor triangular retractors.




PretzelFlex™ 3mm and 5mm

PretzelFlex™ 3mm



Product Name	Product Code	Platform Width	Working Length
PretzelFlex™ 3mm	SHR7301	70mm	231mm
Replacement sterilisation sleeve	SHR7331 (minimum order quantity of 5 units)	(n/a)	304mm

PretzelFlex™ 5mm



Product Name	Product Code	Platform Width	Working Length
PretzelFlex™ 5mm	SHR7001	100mm	285mm
Replacement sterilisation sleeve	SHR7101 (minimum order quantity of 5 units)	(n/a)	365mm

Delivering products designed, developed and manufactured in the UK



 balmer medical

Surgical
Innovations

Balmer Médical SA
Rte de Provence
1426 Concise
info@balmermedical.ch
Tél. +41 24 436 29 07

www.balmermedical.ch

Leaflet Part Number: ALL005 Issue1
Change Request Number: CR3968