


ELECTROSURGERY
COLPOSCOPY



Designed and Manufactured by Prima Medical Ltd

All rights reserved. No part of this publication may be reproduced or transmitted in any form or by any other means, whether electronic or mechanical including by photocopy, recording or information storage and retrieval system, without permission in writing from the publisher. Authorisation to photocopy items for internal or personal use, or the internal or personal use of specific clients is granted by Prima Medical Ltd.

All Products are Latex Free

Copyright © 2006, Prima Medical Ltd.

CONTENTS

Biopsy Punches & Gynae Hook	2
Loop Electrodes (Re-Usable)	3
Fingerswitches & Ball Electrodes (Re-Usable)	4
Loop Electrodes (Single Use)	5
Fingerswitches & Ball Electrodes (Single Use)	6



Prima Medical Ltd

Unit 21, Daniels Industrial Estate, Bath Road, Stroud
Gloucestershire, GL5 3TJ, United Kingdom



+44 (0) 1453 752626






+44 (0) 1453 756750



sales@prima-medical.com



www.prima-medical.com



balmer medical



www.balmermedical.ch

BIOPSY PUNCH



GYNAE HOOK

2



Enlarged for detailed inspection (Not to Scale)

2




(Not to Scale)

BIOPSY PUNCH


1 Biopsy Punch Tishler-Morgan Re-usable pack 1 PMS1000

GYNAE HOOK

2 Gynae Hook 250mm Single Use pack 50 PMS1790



LOOPS




RE-USABLE COLPOSCOPY ELECTRODES

1 Loop Electrode, 8 mm × 8 mm loop, 11cm,	pack 5	PMS1509
2 Loop Electrode, 10 mm × 5 mm loop, 11cm,	pack 5	PMS1510
3 Loop Electrode, 10 mm × 10 mm loop, 11cm,	pack 5	PMS1511
4 Loop Electrode, 12 mm × 8 mm loop, 11cm,	pack 5	PMS1512

5 Loop Electrode, 12 mm × 12 mm loop, 11cm,	pack 5	PMS1501
6 Loop Electrode, 15 mm × 8 mm loop, 11cm,	pack 5	PMS1513
7 Loop Electrode, 15 mm × 15 mm loop, 11cm,	pack 5	PMS1502
8 Loop Electrode, 20 mm × 10 mm loop, 11cm,	pack 5	PMS1514

9 Loop Electrode, 20 mm × 15 mm loop, 11cm,	pack 5	PMS1503
10 Loop Electrode, 25 mm × 10 mm loop, 11cm,	pack 5	PMS1515
11 Loop Electrode, 25 mm × 15 mm loop, 11cm,	pack 5	PMS1517
12 Loop Electrode, 25 mm × 25 mm loop, 11cm,	pack 5	PMS1504

Size may vary by 5%



PRIMA SAFETY FEATURES

Safety guard



Any break in the outer black insulation is highlighted by the yellow insulation layer below

PRIMA SAFETY FEATURES

Safety shield



A slip-on sleeve designed to prevent potential arcing between the electrode and fingerswitch.

PRIMA SAFETY FEATURES

Indicator grooves



Enables sterilisation record to be made on the product. (up to 25 uses)

UNIQUE NUMBER




Enables sterilisation records to be maintained using the appropriate systems.


FINGERSWITCH



BALL ELECTRODES



LOOP ELECTRODES




RE-USABLE FINGERSWITCH

1 100 Plus Fingerswitch, 3m cable terminated with 3 pin plug,
(Autoclavable for a maximum of 100 cycles) pack 5 PMS750

RE-USABLE BALL ELECTRODES

2 3 mm Ball Electrode, 11 cm, pack 3 PMS1408
3 5 mm Ball Electrode, 11 cm, pack 3 PMS1420




LOOPS





SINGLE USE COLPOSCOPY ELECTRODES

Size may vary by 5%

1 Loop Electrode, 10mm x 5mm loop, Pink,	pack 25	PMS1740	5 Loop Electrode, 15mm x 15mm loop, Blue	pack 25	PMS1743	9 Loop Electrode, 20mm x 20mm loop, Beige,	pack 25	PMS1744
2 Loop Electrode, 10mm x 8mm loop, Aqua,	pack 25	PMS1750	6 Loop Electrode, 20mm x 12mm loop, Grey	pack 25	PMS1751	10 Loop Electrode, 5mm x 5mm square loop, Light Green,	pack 25	PMS1735
3 Loop Electrode, 10mm x 10mm loop, Green	pack 25	PMS1742	7 Loop Electrode, 20mm x 13mm loop, Yellow,	pack 25	PMS1752	11 Loop Electrode, 10mm x 4mm square loop, Green	pack 25	PMS1731
4 Loop Electrode, 15mm x 8mm loop, Orange	pack 25	PMS1746	8 Loop Electrode, 20mm x 15mm loop, Red,	pack 25	PMS1745	12 Loop Electrode, 10mm x 8mm square loop, Brown	pack 25	PMS1732






FINGERSWITCH



1

Electrode not included

BALL ELECTRODES




SINGLE USE FINGERSWITCH

1 Fingerswitch with 3 pin plug, 3m cable, boxed 25 PMS255

SINGLE USE BALL ELECTRODES

2 Ball Electrode, 3 mm, 12cm pack 25 PMS1761
3 Ball Electrode, 5 mm, 12cm pack 25 PMS1762



view our online catalogue at

www.prima-medical.com

Call today
to receive your
FREE
Prima catalogue

